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


The image illustrates the steps to add a DAQ Assistant component to a loop in the LabVIEW Express environment. It consists of three main parts:

- Express Panel:** A window titled "Express" showing various functional blocks. A red arrow points to the "Input" block, which is highlighted by a red box containing the text "Wybrać opcję **INPUT**".
- Input Panel:** A window titled "Input" showing a grid of sub-components. A red arrow points to the "DAQ Assist" component, which is highlighted by a red box containing the text "Wybrać funkcję **DAQ Assistant** i przenieść ją do pętli".
- Block Diagram:** A partial view of a LabVIEW block diagram showing a loop structure. A blue arrow points from the "DAQ Assist" component in the "Input" panel to a "DAQ Assistant" component placed inside the loop. A "stop" button is also visible in the diagram.

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Terminal Configuration  
RSE

Timing Settings  
Acquisition Mode: Continuous Samples  
Samples to Read: 100  
Rate (Hz): 1k

$$\text{Odcinek pomiarowy} = \frac{\text{Samples to read}}{\text{Rate (Hz)}} = \frac{100}{1000} = 0,1\text{s}$$

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
The image illustrates the steps to create a graph indicator in LabVIEW. It shows three stages of the process:

- Initial State:** A DAQ Assistant block with a 'data' indicator selected. A red arrow points to the 'data' indicator.
- Context Menu:** A right-click context menu is open over the 'data' indicator. The 'Create' option is highlighted, and its sub-menu is visible, showing 'Graph Indicator' as the selected option. A red arrow points to the 'Create' option, and another red arrow points to 'Graph Indicator'.
- Final State:** The DAQ Assistant block now has a graph indicator connected to the 'data' indicator. A red arrow points to the graph indicator.

Additional elements in the image include a 'stop' button and a 'TF' (True/False) indicator.

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